

Product Datasheet

MEK1 (Thr292) Antibody (orb319572)

Description	Our Anti-MEK1 (Thr292) rabbit polyclonal phosphospecific primary antibody from Biorbyt is produced i...
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic phospho-peptide corresponding to amino acid residues surrounding Thr292 of human MEK1, conjugated to keyhole limpet hemocyanin (KLH).
Target	MEK1 (Thr292)
Preservatives	10 mM HEPES (pH 7.5), 150 mM NaCl, 100 µg per ml rAlbumin and 50% glycerol.
Form/Appearance	Antigen Affinity Purified from Pooled Serum
Storage	Storage at -20°C is recommended, as aliquots may be taken without freeze/thawing due to presence of 50% glycerol. Stable for at least 1 year at -20°C.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	45
Uniprot ID	Q02750
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.