

Product Datasheet

Neurofilament NF-H Antibody (orb319585)

Description	Our Anti-Neurofilament NF-H primary antibody from Biorbyt is chicken polyclonal. It detects all know...
Species/Host	Gallus
Reactivity	Gallus, Mammal
Conjugation	Unconjugated
Tested Applications	ICC, IHC, WB
Immunogen	Purified bovine NF-H
Target	Neurofilament NF-H
Preservatives	PBS + 10 mM NaN ₃ .
Form/Appearance	Total IgY fraction
Storage	Recommended that the undiluted antibody be aliquoted into smaller working volumes (10-30 uL/vial depending on usage) upon arrival and stored long term at -20° C or -80° C, while keeping a working aliquot stored at 4° C for short term. Avoid freeze/thaw cycles. Stable for at least 1 year.
Note	For research use only
Isotype	IgY
Clonality	Polyclonal
MW	200
Uniprot ID	P12036
Dilution Range	WB: 1:50000, IHC: 1:5000-1:25,000, ICC: 1:25,000
Expiration Date	12 months from date of receipt.