

## Product Datasheet

### NMDA NR2A Subunit, N-terminus Antibody (orb319598)

## Description

Our Anti-NMDA NR2A Subunit, N-

antibody from Biorbyt is produced i...

### Species/Host

Rabbit

### Reactivity

Mouse, Rat

### Conjugation

Unconjugated

### Tested

ICC, WB

### Applications

### Immunogen

Synthetic peptide corresponding to amino acid residues from the N-terminal region of the NR2A subunit of rat NMDA receptor, conjugated to keyhole limpet hemocyanin (KLH).

### Target

NMDA NR2A Subunit, N-terminus

### Preservatives

10 mM HEPES (pH 7.5), 150 mM NaCl, 100 µg per ml rAlbumin and 50% glycerol.

### Form/Appearance

Antigen Affinity Purified from Pooled Serum

### Storage

Storage at -20°C is recommended, as aliquots may be taken without freeze/thawing due to presence of 50% glycerol. Stable for at least 1 year at -20°C.

### Note

For research use only

### Isotype

IgG

### Clonality

Polyclonal

### MW

180

### Uniprot ID

[Q00959](#)

### Dilution Range

WB: 1:1000, ICC: 1:500

### Expiration Date

12 months from date of receipt.