

Product Datasheet

Synapsin I Antibody (orb319665)

Description	Our Anti-Synapsin I rabbit polyclonal primary antibody from Biorbyt is produced in-house. It detects...
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IP, WB
Immunogen	Native protein purified from bovine brain.
Target	Synapsin I
Preservatives	10 mM HEPES (pH 7.5), 150 mM NaCl, 100 µg per ml rAlbumin and 50% glycerol.
Form/Appearance	Antigen Affinity Purified from Pooled Serum
Storage	Storage at -20°C is recommended, as aliquots may be taken without freeze/thawing due to presence of 50% glycerol. Stable for at least 1 year at -20°C.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	78
Uniprot ID	P17599
Dilution Range	WB: 1:1000, ICC: 1:1000-1:2000, IP: 1 µl per 200 µg lysate
Expiration Date	12 months from date of receipt.