

Product Datasheet

Synaptojanin 1 Antibody (orb319673)

Description	Our Anti-Synaptojanin 1 primary antibody from Biorbyt is mouse monoclonal. It detects rat Synaptojan...
Species/Host	Mouse
Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Recombinant protein from the C-terminal region of rat synaptojanin 1.
Target	Synaptojanin 1
Preservatives	10 mM HEPES (pH 7.5), 150 mM NaCl, 100 µg per ml rAlbumin and 50% glycerol.
Form/Appearance	Protein G Purified
Storage	Storage at -20°C is recommended, as aliquots may be taken without freeze/thawing due to presence of 50% glycerol. Stable for at least 1 year at -20°C.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	145
Uniprot ID	Q62910
Dilution Range	WB: 1:1000, IHC: 1:100
Expiration Date	12 months from date of receipt.