

Product Datasheet

DNASE1L3 antibody (orb31970)

Description	Peptide affinity purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	DNASE1L3
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	35504
Uniprot ID	Q13609
NCBI	NP_004935.1 , NP_001243489.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

