

Product Datasheet

Ubiquitin C-terminal Hydrolase 1 (UCHL1) Antibody (orb319704)

Description
Our Anti-Ubiquitin C-terminal

antibody from Biorbyt is mouse monoclonal....

Species/Host

Mouse

Reactivity

Bovine, Human, Mouse, Porcine, Rat

Conjugation

Unconjugated

Tested

ICC, IHC, WB

Applications
Immunogen

Recombinant full length human UCHL1 purified from E. coli.

Target

Ubiquitin C-terminal Hydrolase 1 (UCHL1)

Preservatives

PBS + 50% glycerol and 5 mm NaN3

Form/Appearance

Protein G Purified

Storage

Storage at -20°C is recommended, as aliquots may be taken without freeze/thawing due to presence of 50% glycerol. Stable for at least 1 year at -20°C.

Note

For research use only

Isotype

IgG1

Clonality

Monoclonal

MW

24

Uniprot ID

P09936

Dilution Range

WB: 1:5000, IHC: 1:500-1:5000, ICC: 1:500-1:2000

Expiration Date

12 months from date of receipt.