

Product Datasheet

Anti-Vesicular Monoamine Transporter 2, C-terminus (VMAT2) Antibody (orb319710)

Description
Our Anti-Vesicular Monoamine

(VMAT2) sheep polyclonal primary antibody fro...

Species/Host

Sheep

Reactivity

Human, Rat

Conjugation

Unconjugated

Tested

IHC, WB

Applications
Immunogen

Synthetic peptide corresponding to amino acid residues from the intracellular C-terminal region of human VMAT2, conjugated to keyhole limpet hemocyanin (KLH).

Target

Vesicular Monoamine Transporter 2, C-terminus (VMAT2)

Preservatives

10 mM HEPES (pH 7.5), 150 mM NaCl, 100 µg per ml rAlbumin and 50% glycerol.

Form/Appearance

Antigen Affinity Purified from Pooled Serum

Storage

Storage at -20°C is recommended, as aliquots may be taken without freeze/thawing due to presence of 50% glycerol. Stable for at least 1 year at -20°C.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

MW

57

Uniprot ID

Q05940

Dilution Range

WB: 1:1000, IHC: 1:200

Expiration Date

12 months from date of receipt.