

Product Datasheet

ALDH1L1 Antibody (orb319716)

Description	ALDH1L1 Antibody
Species/Host	Mouse
Reactivity	Bovine, Gallus, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, WB
Immunogen	Amino acids 1- 401 of human ALDH1L1
Form/Appearance	Affinity Purified
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Isotype	Mouse IgG
Uniprot ID	075891
Entrez	10840
Expiration Date	12 months from date of receipt.