

Product Datasheet

Glutamine Synthetase Antibody (orb319768)

Description	Glutamine Synthetase Antibody
Species/Host	Rabbit
Reactivity	Bovine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Recombinant Protein
Form/Appearance	Whole serum
Storage	Aliquot and store at -20°C. Avoid freeze-thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Uniprot ID	P15104
Entrez	2752
Expiration Date	12 months from date of receipt.