



Product Datasheet

NFATC1 antibody (orb31985)

Description	Peptide affinity purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Synthetic Peptide
Target	NFATC1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	101243
Uniprot ID	O95644
NCBI	NP_001265601.1 , NP_001265602.1 , NP_001265604.1 , NP_765975.1 , NP_006153.2 , NP_765976.1 , NP_001265599.1 , NP_001265598.1 , NP_765978.1 , NP_765977.1
Dilution Range	WB: 1:1000, IHC-P: 1:10-50
Expiration Date	12 months from date of receipt.



Western blot analysis of NCI-H292 cell l...