

Product Datasheet

Simeprevir(TMC-435) (orb320240)

Description	Simeprevir (TMC 435350) is a potent inhibitor of HCV NS3/4A protease (K _i = 0.36 nM) and viral replic...
Note	For research use only
Formula	C ₃₈ H ₄₇ N ₅ O ₇ S ₂
MW	749.94
CAS Number	[923604-59-5]
SMILES	<chem>COC1=CC=C2C(=CC(C3SC=C(C(C)C)N=3)=NC2=C1C)O[C@H]1C[C@H]2C(N(CCCCC=C[C@@H]3C[C@]3(NC(=O)[C@@H]2C1)C(NS(C1CC1)(=O)=O)=O)C)=O</chem> t:34,&1:22,24,33,35,39
Expiration Date	12 months from date of receipt.