

Product Datasheet

C4orf18 antibody (orb32051)

Description	Peptide affinity purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	GASK1B
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	57552
Uniprot ID	Q6UWH4
NCBI	NP_001121896.1 , NP_001026870.2 , NP_057697.2
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

