

Product Datasheet

LAMC2 antibody (orb32069)

Description	Peptide affinity purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	LAMC2
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	130976
Uniprot ID	Q13753
NCBI	NP_005553.2 , NP_061486.2
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

