



Product Datasheet

MYB antibody (orb32075)

Description	Peptide affinity purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Synthetic Peptide
Target	MYB
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	72341
Uniprot ID	P10242
NCBI	NP_001155130.1 , NP_001155128.1 , NP_001155132.1 , NP_001155129.1 , NP_005366.2 , NP_001123645.1 , NP_001155131.1 , NP_001123644.1
Dilution Range	WB: 1:1000, IHC-P: 1:10-50
Expiration Date	12 months from date of receipt.



Western blot analysis of NCI-H292 cell l...