

Product Datasheet

Laminin Peptide YIGSR-NH2 (orb320879)

Description	Laminin Peptide YIGSR-NH2
Storage	Stability: Store in a cool, dry place.. Storage Temp: -20°C
Note	For research use only
Formula	C ₂₆ H ₄₃ N ₉ O ₇
Purity	≥98%
MW	593.67
SMILES	<chem>[H]N[C@@H](CC1=CC=C(O)C=C1)C(=O)N[C@@H]([C@@H](C)CC)C(=O)NCC(=O)N[C@@H](CO)C(=O)N[C@@H](CCCNC(=N)NN)C=O</chem>
Solubility (25°C)	Soluble in water.
Hazard Information	UN Number: Not dangerous goods
Expiration Date	12 months from date of receipt.