

Product Datasheet

Galanin (orb321070)

Description	Galanin
Storage	Storage Temp: -20°C
Note	For research use only
Formula	C ₁₃₉ H ₂₁₀ N ₄₂ O ₄₃
Purity	≥98%
MW	3157.44
CAS Number	119418-04-1
SMILES	<chem>CC(C)CC(C(=O)NC(C(C)O)C(=O)NC(CO)C(=O)O)NC(=O)CNC(=O)C(CC(=O)N)NC(=O)C(CCCCN)NC(=O)C(CC(=O)O)NC(=O)C(CO)NC(=O)C(CC1=CC=CC=C1)NC(=O)C(CO)NC(=O)C(CCCN(C=N)N)NC(=O)C(CC2=CNC=N2)NC(=O)C(CC(=O)N)NC(=O)CNC(=O)C(C(C)C)NC(=O)C(C)NC(=O)C(C3=CNC=N3)NC(=O)C4CCCN4C(=O)CNC(=O)C(CC(C)C)NC(=O)C(CC(C)C)NC(=O)C(CC5=CC=C(C=C5)O)NC(=O)CNC(=O)C(C)NC(=O)C(CO)NC(=O)C(CC(=O)N)NC(=O)C(CC(C)C)NC(=O)C(C(C)O)NC(=O)C(CC6=CNC7=CC=CC=C76)NC(=O)CN</chem>
Solubility (25°C)	Soluble in 5% acetic acid.
Hazard Information	UN Number: Not dangerous goods
Expiration Date	12 months from date of receipt.