

Product Datasheet

PPP2R5C antibody (orb32122)

Description	Peptide affinity purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	PPP2R5C
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	61061
Uniprot ID	Q13362
NCBI	NP_848702.1 , NP_001155198.1 , NP_848701.1 , NP_001155197.1 , NP_002710.2
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

