

Product Datasheet

Ac-IEAR-pNA (orb321368)

Description	Ac-IEAR-pNA
Storage	Storage Temp: -20°C
Note	For research use only
Formula	C ₂₈ H ₄₃ N ₉ O ₉
Purity	≥95%
MW	649.7
SMILES	<chem>CC[C@@H](C)[C@@H](C(=O)N[C@@H](CCC(=O)O)C(=O)N[C@@H](C)C(=O)N[C@@H](CCCNC(=N)N)C(=O)Nc1ccc(cc1)[N+](=O)[O-])NC(=O)C</chem>
Solubility (25°C)	Soluble in water.
Expiration Date	12 months from date of receipt.