

Product Datasheet

Adipokinetic Hormone (orb321382)

Description	Adipokinetic Hormone
Storage	Storage Temp: -20°C
Note	For research use only
Formula	C ₄₄ H ₆₀ N ₁₀ O ₁₂
Purity	≥95%
MW	921.03
SMILES	<chem>CC(C)C[C@H](NC(=O)[C@@H]1CCC(=O)N1)C(=O)N[C@@H]([C@@H](C)O)C(=O)N[C@@H](CC1=CC=CC=C1)C(=O)N[C@@H]([C@@H](C)O)C(=O)N[C@@H](CO)C(=O)N[C@@H](CC1=CNC2=C1C=CC=C2)C(=O)NCC(N)=O</chem>
Solubility (25°C)	Soluble in water.
Hazard Information	UN Number: Not dangerous goods
Expiration Date	12 months from date of receipt.