

Product Datasheet

Angiotensin Converting Enzyme Inhibitor Peptide (orb321411)

Description	Angiotensin Converting Enzyme Inhibitor Peptide
Storage	Storage Temp: -20°C
Note	For research use only
Application notes	White to off white powder
Formula	C ₅₃ H ₇₇ N ₁₄ O ₁₂
Purity	≥95%
MW	1102.29
SMILES	<chem>[H]N1[C@@H](CCC1=O)C=[O]/N[C@@H](CC1=CNC2=C1C=CC=C2)C(=O)N1CCC[C@H]1C(=O)N[C@@H](CCCNC(N)=N)C(=O)N1CCC[C@H]1C(=O)N[C@@H](CCC(N)=O)C(=O)N[C@@H]([C@@H](C)CC)C(=O)N1CCC[C@H]1C(=O)N1CCC[C@H]1C(O)=O</chem>
Solubility (25°C)	Soluble in water.
Hazard Information	UN Number: Not dangerous goods
Expiration Date	12 months from date of receipt.