

Product Datasheet

[Des-Asp1]-Angiotensin I, human (orb321412)

Description	[Des-Asp1]-Angiotensin I, human
Storage	Storage Temp: -20°C
Note	For research use only
Application notes	White to off white powder
Formula	C ₅₈ H ₈₄ N ₁₆ O ₁₁
Purity	≥95%
MW	1181.42
SMILES	<chem>[H]N[C@@H](CCCNC(N)=N)C(=O)N[C@@H](C(C)C)C(=O)N[C@@H](CC1=CC=C(O)C=C1)C(=O)N[C@@H]([C@@H](C)CC)C(=O)N[C@@H](CC1=CNC=N1)C(=O)N1CCC[C@H]1C(=O)N[C@@H](CC1=CC=CC=C1)C(=O)N[C@@H](CC1=CNC=N1)C(=O)N[C@@H](CC(C)C)C(=O)O</chem>
Solubility (25°C)	Soluble in water.
Expiration Date	12 months from date of receipt.