

Product Datasheet

Bradykinin Potentiator B (orb321458)

Description	Bradykinin Potentiator B
Storage	Storage Temp: -20°C
Note	For research use only
Formula	C ₅₆ H ₉₁ N ₁₅ O ₁₃
Purity	≥95%
MW	1182.46
SMILES	<chem>O=C(N1[C@H](C(=O)O)CCC1)[C@H]6N(C(=O)[C@@H](NC(=O)[C@@H](NC(=O)[C@H]5N(C(=O)[C@@H](NC(=O)[C@H]4N(C(=O)[C@H]3N(C(=O)[C@@H](NC(=O)CNC(=O)[C@H]2NC(=O)CC2)CC(C)C)CCC3)CCC4)CCCN=C(\N)N)CCC5)CCCCN)[C@@H](C)CC)CCC6</chem>
Solubility (25°C)	Soluble in water.
Hazard Information	UN Number: Not dangerous goods
Expiration Date	12 months from date of receipt.