

Product Datasheet

C-Telopeptide (orb321505)

Description	C-Telopeptide
Storage	Storage Temp: -20°C
Note	For research use only
Formula	C ₃₄ H ₅₆ N ₁₄ O ₁₃
Purity	≥95%
MW	868.91
SMILES	<chem>O=C(O)[C@@H](NC(=O)CNC(=O)CNC(=O)[C@@H](NC(=O)[C@@H](NC(=O)[C@@H](NC(=O)[C@@H](N)CCC(=O)O)CCCN)C)Cc1cncn1)CC(=O)O)CCC/N=C(\N)N</chem>
Solubility (25°C)	Soluble in water.
Hazard Information	UN Number: Not dangerous goods
Expiration Date	12 months from date of receipt.