



Product Datasheet

SEC22B antibody (orb32155)

Description	Rabbit polyclonal antibody to SEC22B functions as a v-SNARE of transport vesicles and plays a role i...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Synthetic Peptide
Target	SEC22B
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	24741
Uniprot ID	O75396
NCBI	NP_004883.3
Dilution Range	WB: 1:1000, IHC-P: 1:10-50
Expiration Date	12 months from date of receipt.



Western blot analysis of SEC22B antibody...