

Product Datasheet

G250.A2 Peptide (orb321555)

Description	G250.A2 Peptide
Storage	Storage Temp: -20°C
Note	For research use only
Formula	$C_{45}H_{72}N_{14}O_{12}$
Purity	≥95%
MW	1001.16
SMILES	<chem>[H]N[C@@H](CC1=CNC=N1)C(=O)N[C@@H](CC(C)C)C(=O)N[C@@H](CO)C(=O)N[C@@H]([C@@H](CO)C(=O)N[C@@H](C)C(=O)N[C@@H](CC1=CC=CC=C1)C(=O)N[C@@H](C)C(=O)N[C@@H](CCCNC(N)=N)C(=O)N[C@@H](C(C)C)C(O)=O</chem>
Expiration Date	12 months from date of receipt.