

Product Datasheet

Peptide 401 (orb321632)

Description	Peptide 401
Storage	Storage Temp: -20°C
Note	For research use only
Formula	C ₁₁₀ H ₁₈₈ N ₄₀ O ₂₄ S ₄
Purity	≥95%
MW	2583.28
CAS Number	32908-73-9
SMILES	<chem>CCC(C)C1C(=O)NC(CSSCC(C(=O)NC(C(=O)NC2CSSCC(C(=O)NC(C(=O)NC(C(=O)N1)CCCCN)CCCNC(=N)N)NC(=O)C(NC(=O)C(NC(=O)C3CCCCN3C(=O)C(NC(=O)C(NC(=O)C(NC(=O)C(NC(=O)C(NC2=O)CCCCN)CCCNC(=N)N)CC4C=NC=N4)C(C)C)C(C)CC)CCCCN)CC5C=NC=N5)C(C)C)CC(=O)N)NC(=O)C(CCCCN)NC(=O)C(C(C)CC)N)C(=O)NCC(=O)NC(CC)CCN)C(=O)NC(CC(=O)N)C(=O)N</chem>
Solubility (25°C)	Soluble in water.
Expiration Date	12 months from date of receipt.