

Product Datasheet

P1 Peptide (orb321666)

Description	P1 Peptide
Storage	Storage Temp: -20°C
Note	For research use only
Formula	C ₃₉ H ₆₃ N ₁₁ O ₁₃
Purity	≥95%
MW	894
SMILES	<chem>[H]NCC(=O)N[C@@H](CO)C(=O)N[C@@H](CC1=CC=CC=C1)C(=O)N[C@@H](CC(C)C)C(=O)N[C@@H](C(C)C)C(=O)N[C@@H](CCNC(N)=N)C(=O)N[C@@H](CCC(O)=O)C(=O)N[C@@H](CO)C(O)=O</chem>
Expiration Date	12 months from date of receipt.