

Product Datasheet

P55-TNFR Peptide (orb321667)

Description	P55-TNFR Peptide
Storage	Storage Temp: -20°C
Note	For research use only
Formula	C ₅₇ H ₉₅ N ₁₅ O ₂₂
Purity	≥95%
MW	1342.48
SMILES	[H]N[C@@H](CC(C)C)C(=O)N1CCC[C@H]1C(=O)N[C@@H](CCC(N)=O)C(=O)N[C@@H]([C@@H](C)CC)C(=O)N[C@@H](CCC(O)=O)C(=O)N[C@@H](CC(N)=O)C(=O)N[C@@H](C(C)C)C(=O)N[C@@H](CCCCN)C(=O)NCC(=O)N[C@@H]([C@@H](C)O)C(=O)N[C@@H](CCC(O)=O)C(=O)N[C@@H](CC(O)=O)C(O)=O
Solubility (25°C)	Soluble in water.
Expiration Date	12 months from date of receipt.