

Product Datasheet

Phosphate Acceptor Peptide (orb321690)

Description	Phosphate Acceptor Peptide
Storage	Storage Temp: -20°C
Note	For research use only
Application notes	White to off white powder
Formula	$C_{41}H_{74}N_{16}O_{11}$
Purity	≥95%
MW	967.15
SMILES	<chem>[H]N[C@@H](CCCNC(N)=N)C(=O)N[C@@H](CCCNC(N)=N)C(=O)N[C@@H](CCCCN)C(=O)N[C@@H](C)C(=O)N[C@@H](CO)C(=O)NCC(=O)N1CC[C@H]1C(=O)N1CCC[C@H]1C(=O)N[C@@H](C(C)C)C(O)=O</chem>
Solubility (25°C)	Soluble in water.
Hazard Information	UN Number: Not dangerous goods
Expiration Date	12 months from date of receipt.