

Product Datasheet

Renin Inhibitor Peptide (orb321703)

Description	Renin Inhibitor Peptide
Storage	Storage Temp: -20°C
Note	For research use only
Formula	C ₅₂ H ₇₃ N ₁₃ O ₉
Purity	≥95%
MW	1024.24
SMILES	<chem>CC(C)C[C@H](NC(=O)[C@H](CC1=CNC=N1)NC(=O)[C@H](CC1=CC=CC=C1)NC(=O)[C@@H]1CCCN1C(=O)[C@@H](N)CC1=CNC=N1)C(=O)N[C@H](CC(C)C)C(=O)N[C@@H](C(C)C)C(=O)N[C@@H](CC1=CC=C(N)C=C1)C(O)=O</chem>
Solubility (25°C)	Soluble in water.
Expiration Date	12 months from date of receipt.