

Terlipressin Acetate (orb321756)

Description	Terlipressin Acetate
Storage	Storage Temp: -20°C
Note	For research use only
Application notes	White to off white powder
Formula	C ₅₂ H ₇₄ N ₁₆ O ₁₅ S ₂
Purity	≥95%
MW	1227.37
CAS Number	14636-12-5
SMILES	C1CC(N(C1)C(=O)C2CSSCC(C(=O)NC(C(=O)NC(C(=O)NC(C(=O)NC(C(=O)N2)CC(=O)N)CCC(=O)N)CC3=CC=CC=C3)CC4=CC=C(C=C4)O)NC(=O)CNC(=O)CNC(=O)CN)C(=O)NC(CCCCN)C(=O)NCC(=O)N
Solubility (25°C)	Insoluble in water. Soluble in DMSO.
Hazard Information	UN Number: Not dangerous goods
Expiration Date	12 months from date of receipt.