

Product Datasheet

Thrombin Receptor Agonist Peptide (orb321758)

Description	Thrombin Receptor Agonist Peptide
Storage	Storage Temp: -20°C
Note	For research use only
Application notes	White to off white powder
Formula	C ₈₁ H ₁₁₈ N ₂₀ O ₂₃
Purity	≥95%
MW	1739.96
CAS Number	137339-65-2
SMILES	<chem>CC(C)CC(C(=O)NC(CC(C)C)C(=O)NC(CCCNC(=N)N)C(=O)NC(CC(=O)N)C(=O)N1CCCC1C(=O)NC(CC(=O)N)C(=O)NC(CC(=O)O)C(=O)NC(CCCCN)C(=O)NC(CC2=CC=C(C=C2)O)C(=O)NC(CCC(=O)O)C(=O)N3CCCC3C(=O)NC(CC4=CC=CC=C4)C(=O)O)NC(=O)C(CC5=CC=CC=C5)NC(=O)C(CO)N</chem>
Solubility (25°C)	Soluble in water (0.9 mg/mL).
Hazard Information	UN Number: Not dangerous goods
Expiration Date	12 months from date of receipt.