

Product Datasheet

Universal Tetanus Toxin Epitope P2 (830-844) (orb321766)

Description	Universal Tetanus Toxin Epitope P2 (830-844)
Storage	Storage Temp: -20°C
Note	For research use only
Formula	$C_{80}H_{129}N_{19}O_{23}$
Purity	≥95%
MW	1725.03
SMILES	<chem>[H]N[C@@H](CCC(N)=O)C(=O)N[C@@H](CC1=CC=C(O)C=C1)C(=O)N[C@@H]([C@@H](C)CC)C(=O)N[C@@H](CCCCN)C(=O)N[C@@H](C)C(=O)N[C@@H](CC(N)=O)C(=O)N[C@@H](CO)C(=O)N[C@@H](CCCCN)C(=O)N[C@@H](CC1=CC=CC=C1)C(=O)N[C@@H]([C@@H](C)CC)C(=O)NC(=O)N[C@@H]([C@@H](C)CC)C(=O)N[C@@H]([C@@H](C)O)C(=O)N[C@@H](CCC(O)=O)C(=O)N[C@@H](CC(C)C)C(=O)O</chem>
Hazard Information	UN Number: Not dangerous goods
Expiration Date	12 months from date of receipt.