



Product Datasheet

SOX9 antibody (orb32198)

Description	Rabbit polyclonal antibody to UTF1 is expressed in cells of the inner cell mass, and the epiblast t...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	Synthetic Peptide
Target	SOX9
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	56137
Uniprot ID	P48436
NCBI	NP_000337.1
Dilution Range	WB: 1:1000, FACS: 1:10-50
Expiration Date	12 months from date of receipt.



Flow
cytometric
analysis of
HepG2 cells
...