



Product Datasheet

TIMM23 antibody (orb32202)

Description	Peptide affinity purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	TIMM23
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	21943
Uniprot ID	O14925
NCBI	NP_006318.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.



WB analysis
of MCF-7
A431 cell line
(fro...