


Product Datasheet

Tocotrienol (orb322107)

Description	Tocotrienol
Storage	Storage Temp: -20°C
Note	For research use only
Application notes	yellow liquid
Formula	C ₂₉ H ₄₄ O ₂
Purity	≥98%
MW	424.66
CAS Number	58864-81-6
SMILES	<chem>CC(=CCC/C(=C/CC/C(=C/CC[C@]1(C)CCc2c(C)c(c(C)c(C)c2O1)O)/C)/C)C</chem>
Solubility (25°C)	Soluble in ethanol, hexanes, DMSO, DMF
Hazard Information	UN Number: Not dangerous goods
Expiration Date	12 months from date of receipt.



Chemical
structure of
Tocotrienol