


Product Datasheet

PNKP antibody (orb32218)

Description	Peptide affinity purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	PNKP
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Concentration	0.5 mg/ml
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	57076
Uniprot ID	Q96T60
NCBI	NP_009185.2
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of 293 cell lysate...

