


Product Datasheet

TM3-0089 Sodium Monohydrate (orb322295)

Description	TM3-0089 Sodium Monohydrate
Storage	Storage Temp: -20°C
Note	For research use only
Application notes	White to off white powder
Formula	$C_{21}H_{20}FN_2O_4S Na H_2O$
Purity	≥98%
MW	456.46
CAS Number	844639-57-2
SMILES	[Na+].[H]O[H].CN(C1CCC2C(C1)C1=CC=CC=C1N2CC([O-])=O)S(=O)(=O)C1=CC=C(F)C=C1
Hazard Information	UN Number: Not dangerous goods
Expiration Date	12 months from date of receipt.



Chemical
structure of
TM3-0089
Sodium Mo...