

Product Datasheet

PF14_0076 antibody (orb32250)

Description	PF14_0076 antibody
Species/Host	Rabbit
Reactivity	Bacteria
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant Plasmodium falciparum Plasmepsin-1 protein (125-452AA)
Target	PF14_0076
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q7KQM4
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.