


Product Datasheet

Prostaglandin E1 (orb322660)

Description	Prostaglandin E1
Storage	Storage Temp: -20°C
Note	For research use only
Formula	$C_{20}H_{34}O_5$
Purity	≥98%
MW	354.48
CAS Number	745-65-3
SMILES	<chem>CCCCC(C=CC1C(CC(=O)C1CCCCCCC(=O)O)O)O</chem>
Solubility (25°C)	Soluble in ethanol (5mg/mL) or acetone (10mg/mL).
Hazard Information	Non-Toxic
Expiration Date	12 months from date of receipt.



Chemical
structure of
Prostaglandin
E1