


Product Datasheet

Aliskiren Hemifumarate (orb322741)

Description	Aliskiren Hemifumarate
Storage	Storage Temp: -20°C
Note	For research use only
Formula	$C_{30}H_{53}N_3O_6 \cdot 1/2 C_4H_4O_4$
Purity	≥98%
MW	1219.59
CAS Number	173334-58-2
SMILES	<chem>CC(C)[C@@H](Cc1ccc(c(c1)OCCOC)OC)C[C@@H]([C@H](C[C@@H](C(C)C)C(=O)NCC(C)(C)C(=O)N)O)N.CC(C)[C@@H](Cc1ccc(c(c1)OCCOC)OC)C[C@@H]([C@H](C[C@@H](C(C)C)C(=O)NCC(C)(C)C(=O)N)O)N.C(=O)C(=O)O</chem>
Solubility (25°C)	Soluble in water (≥122 mg/mL), ethanol ≥122 mg/mL, and DMSO (≥122 mg/mL).
Hazard Information	UN Number: Not dangerous goods
Expiration Date	12 months from date of receipt.



Chemical
structure of
Aliskiren
Hemifuma...