

Product Datasheet

Collagen III antibody (orb322981)

Description

Collagen III antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, ICC, IF, IHC-P

Applications

Immunogen

KLH conjugated synthetic peptide derived from human Collagen III. Please contact us for the exact immunogen sequence. The peptide is available as orb374940.

Target

Collagen III

Form/Appearance

10 mM PBS, 0.02% sodium azide

Concentration

- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml

Storage

Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

Purity

Polyclonal antibodies are purified by peptide affinity chromatography

MW

117 kDa

Uniprot ID

[P08121](#), [P13941](#), [P02461](#)

Entrez



[1281](#)

Dilution Range


IHC-P: 1:100-500, IF/ICC: 1:100-500 (based on 0.5mg/ml)

Expiration Date


12 months from date of receipt.



Immunofluorescence analysis of rat stoma...



IHC-P image of rat spinal cord tissue us...



Immunohistochemical staining of paraffin...