

Product Datasheet

Ly-6G antibody (orb322983)

Description

Ly-6G antibody

Species/Host

Rabbit

Conjugation

Unconjugated

Tested
Applications

ELISA, ICC, IF, IHC-P

Immunogen

KLH conjugated synthetic peptide derived from human Ly-6G. Please contact us for the exact immunogen sequence. The peptide is available as orb374942.

Target

Ly-6G

Form/Appearance

10 mM PBS, 0.02% sodium azide

Concentration

- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml

Storage

Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

Purity

Polyclonal antibodies are purified by peptide affinity chromatography

MW



14 kDa

Uniprot ID
P35461
NCBI
081642, 027366, 11
Entrez
546644
Dilution Range


IHC-P: 1:100-500, IF/ICC: 1:100-500 (based on 0.5mg/ml)

Expiration Date


12 months from date of receipt.



Immunofluorescence analysis of mouse tissue...



Immunohistochemical staining of paraffin-embedded tissue...



IHC-P staining of human endometrial cancer tissue...