



Product Datasheet

IL8 antibody (orb322984)


Description	IL8 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FACS, IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human IL8. Please contact us for the exact immunogen sequence. The peptide is available as orb374943.
Target	IL8
Form/Appearance	10 mM PBS, 0.02% sodium azide
Concentration	- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Polyclonal antibodies are purified by peptide affinity chromatography
Uniprot ID	P10145
Entrez	3576
Dilution Range	IHC-P: 1:100-800, IF/ICC: 1:100-800 (based on 0.5mg/ml)
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of mouse brain tissue showing IL8 expression.



Immunofluorescence image of rat brain tissue showing IL8 expression.



IHC-P image of mouse brain tissue showing IL8 expression using the IL8 antibody (orb322984).