


Product Datasheet

ICOS ligand antibody (orb323034)

Description	Goat polyclonal antibody to ICOS ligand
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	ICOS ligand / B7-H2
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 33kDa band observed in Human Lymph Node lysates (calculated MW of 33.3kDa according to NP_056074.1). Recommended concentration: 1-3µg/ml.
Protein Sequence	SQTGNDIGERDKIT
Clonality	Polyclonal
MW	33.3; 34.2; 20.1; 23.7
Entrez	23308
Dilution Range	ELISA: 1:8000, WB: 1-3 µg/ml
Expiration Date	12 months from date of receipt.



orb323034
(1µg/ml)
staining of
Human Ly...