

Product Datasheet

CO1A2 antibody (orb32306)

Description	Rabbit polyclonal antibody to CO1A2
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	COL2A1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	141785
Uniprot ID	P02458
NCBI	NP_149162.2, NP_001835.3
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

