



Product Datasheet

Acetyl Lysine antibody (orb323086)

Description	Mouse monoclonal antibody to Acetyl Lysine
Species/Host	Mouse
Reactivity	All
Conjugation	Unconjugated
Tested Applications	IH, WB
Immunogen	Purified protein corresponding to Acetyl Lysine.
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Monoclonal
Source	Mouse
Dilution Range	WB: 1:1000-1:2000, IHC-P: 1:200-1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of Acetyl Lysine e...



Immunohistochemical analysis of Acetyl L...