



Product Datasheet

Alpha-tubulin (AcK40) antibody (orb323091)

Description	Mouse monoclonal antibody to Alpha-tubulin (Acetyl-K40)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IH, WB
Immunogen	KLH-conjugated synthetic acetylated peptide corresponding to residues surrounding K40 of human Alpha-tubulin protein. The exact sequence is proprietary.
Target	TUBA1B
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Monoclonal
Source	Mouse
Uniprot ID	P68363
Entrez	7277
Dilution Range	WB: 1:1000-1:2000, IHC-P: 1:100-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of Alpha-tubulin (...)



Immunohistochemical analysis of Alpha-tu...