



Product Datasheet

Tubulin beta 1 antibody (orb323098)

Description	Mouse monoclonal antibody to Tubulin beta 1
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IH, WB
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence of human Beta1-tubulin. The exact sequence is proprietary.
Target	TUBB1
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Monoclonal
Source	Mouse
Uniprot ID	Q9H4B7 , A2AQ07
Entrez	545486 , 81027
Dilution Range	WB: 1:5000-1:10000, IHC-P: 1:200-1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of Beta1-tubulin e...



Immunohistochemical analysis of Beta1-tu...